

Tork Mini Jumbo Bath Tissue Dispenser, Stainless Steel

465500



Description

Tork Image Design Mini Jumbo Bath Tissue Roll Dispenser with Reserve system features a gently brushed stainless steel with fine lines and subtle curves provide clean lines and smooth surfaces. Roll break prevents paper over spin. Stub roll feature allows for full roll consumption. Large capacity rolls allow service for high-traffic areas without adding extra dispensers. Roll is covered to minimize cross contamination.

- Smart, stub-roll holder allows for full roll consumption to minimize waste
- Optimal roll brake avoids free spinning rolls and therefore reduces consumption
- Hygienic: only touch the sheet being used
- ADA compliant
- Tork Vision Cleaning Enabled

Product Certifications



ADA Compliant



Product Details

Product height (without packaging)	10 in
Dispenser Material	Metal/Plastic
Product width (without packaging)	14 in
Product length/depth (without packaging)	5 in
System	T2
Color	Stainless, Stainless steel look

Shipping Data

	Consumer Units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	10073286631258	10073286631258	7322540845358
Packaging Material	Box	Carton	-
Pieces	1	1 (1 CON)	160 (160 TRP)
Height	6.26 in	6.26 in	100.16 in
Length	10.79 in	14.76 in	44.29 in
Width	14.76 in	10.79 in	36.34 in
Gross Weight	3.09 lb	3.71 lb	592.96 lb
Net Weight	3.09 lb	3.09 lb	493.77 lb
Volume	0.58 ft	0.58 ft	93.29 ft
Layers Per Pallet	-	-	16
TRP Per Layer	-	-	10



Tork Mini Jumbo Bath Tissue Dispenser, Stainless Steel

465500

Compatible Products



ADVC JUMBO BATH 2-PLY 12/751
11020602



ADVC JUMBO BATH 1 PLY 12/1200
12013903



ADVC JUMBO BATH 2-PLY 12/751
12024402