



Description

The Tork Skincare Dispenser in Elevation Design fits any location and can be used with a wide range of Tork soaps and sanitizers in sealed bottles. It is Easy to Use certified and provides good hand hygiene for all users. Elevation dispensers have a functional and modern design, that makes a lasting impression on your guests.

- ADA compliant
- Proven Easy to use design by third party certificate -promotes good hand hygiene for all users, even children and elderly
- Reduce time for cleaning staff: Certified effortless refilling ¹ in less than 10 sec ²
- High-capacity reducing waste: Helps reduce consumption by up to 50% ³
- Helps secure good hand hygiene:Factory-sealed bottle with a new pump for every refill helps reduce risk of contamination and protect the formulations all the way to the hands.
- Carbon neutral – reduced and offset
- Carbon neutral – produced with purchased certified renewable electricity and remaining carbon emissions are offset with credits from climate projects

Product Certifications



ADA Compliant



ClimatePartner
certified product
climate-id.com/9VIUDN

Product Details

Product height (without packaging)	11.5 in
Dispenser Material	plastic
Product width (without packaging)	4.45 in
Product length/depth (without packaging)	4.49 in
System	S4
Color	White, White

Shipping Data

	Consumer Units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	73286627322	10073286627329	7322540845617
Packaging Material	Box	Carton	-
Pieces	1	4 (4 CON)	448 (112 TRP)
Height	11.93 in	13 in	91 in
Length	4.65 in	10.13 in	40.5 in
Width	4.76 in	10.13 in	40.5 in
Gross Weight		4.84 lb	542.48 lb
Net Weight	0.8 lb	3.2 lb	358.41 lb
Volume	0.15 ft	0.77 ft	86.38 ft
Layers Per Pallet	-	-	7
TRP Per Layer	-	-	16

1. Certified by the Swedish Rheumatism association, tested according to The Design for All Test, published in the Swedish Design Research Journal, 2011.
2. Essity test: Using Tork Foam Soap 0,4 ml/push vs Tork Liquid Soap 1,0ml/push in 1000ml refill in an Elevation dispenser

Compatible Products



TORK ODOR CONTROL SOAP 6/1L
400020



PREM ALCOHOL FOAM HAND SANI 12/18OZ
400022



PINK LUXURY SOAP 6/1L
400032

PREM ALCOHOL FOAM HAND SANI 6/1L
400217