



Think ahead.

Tork Skincare Dispenser With Intuition® Sensor

571608



Description

Tork Automatic Soap and Sanitizer Dispenser in Elevation Design allows you to make a good impression on your guests and can be used with a wide range of Tork skincare products in sealed bottles. Its touch free dispensing secures good hand hygiene for all users. Elevation dispensers have a functional and modern design that makes a lasting impression on your guests.

- Tork Vision Cleaning Enabled
- Sensor driven, touch-free dispensing reduces the risk of cross contamination
- Reduce time for cleaning staff: Certified effortless refilling¹ in less than 10 sec²
- High-capacity reducing waste: Helps reduce consumption by up to 50%³
- Helps secure good hand hygiene: Factory-sealed bottle with a new pump for every refill helps reduce risk of contamination and protect the formulations all the way to the hands.
- Reliable performance and intuitive maintenance with an LED refill indicator that shows when maintenance is needed, delivers peace of mind.
- Sensor
- Elevation
- Flexibility

Product Certifications



Elevation
Design



Product Details

Product height (without packaging)	10.9 in
Dispenser Material	plastic
Product width (without packaging)	4.41 in
Product length/depth (without packaging)	5.04 in
System	S4
Color	Black

Shipping Data

	Consumer Units (CON)	Transport unit (TRP)	Pallet (PAL)
EAN	73286627407	10073286627404	7322540845648
Packaging Material	Box	Carton	-
Pieces	1	4 (4 CON)	288 (72 TRP)
Height	5.51 in	13.5 in	81 in
Length	6.1 in	12.52 in	50 in
Width	11.81 in	11.61 in	34.88 in
Gross Weight	1.5 lb	8.7 lb	626.36 lb
Net Weight	1.5 lb	6 lb	432.01 lb
Volume	0.23 ft	1.14 ft	81.74 ft
Layers Per Pallet	-	-	6
TRP Per Layer	-	-	12

1. Essity test: Using Tork Foam Soap 0,4 ml/push vs Tork Liquid Soap 1,0ml/push in 1000ml refill in an Elevation dispenser



Think ahead.

Tork Skincare Dispenser With Intuition® Sensor

571608

Compatible Products



TORK ODOR CONTROL SOAP 6/1L
400020



PINK LUXURY SOAP 6/1L
400032



PREM ALCOHOL FOAM HAND SANI 6/1L
400217

CLARITY HAND WASHING FOAM SOAP 6/1L
401800